



Chapter 8

16S rRNA Gene Analysis with QIIME2

Michael Hall and Robert G. Beiko

Abstract

Microbial marker-gene sequence data can be used to generate comprehensive taxonomic profiles of the microorganisms present in a given community and for other community diversity analyses. The process of going from raw gene sequences to taxonomic profiles or diversity measures involves a series of data transformations performed by numerous computational tools. This includes tools for sequence quality checking, denoising, taxonomic classification, alignment, and phylogenetic tree building. In this chapter, we demonstrate how the Quantitative Insights Into Microbial Ecology version 2 (QIIME2) software suite can simplify 16S rRNA marker-gene analysis. We walk through an example data set extracted from the guts of bumblebees in order to show how QIIME2 can transform raw sequences into taxonomic bar plots, phylogenetic trees, principal co-ordinates analyses, and other visualizations of microbial diversity.

Key words Microbial ecology, Marker gene, 16S rRNA gene, QIIME, Bioinformatics

1 Introduction

Molecular techniques give us the ability to characterize microorganisms and gain insight into the important biological processes that they drive. Modern high-throughput methods allow for the interrogation of entire communities of microorganisms in parallel. One such method is marker gene analysis. In this type of microbial community analysis, a phylogenetically informative and universal gene (or gene fragment) is isolated and amplified with the polymerase chain reaction (PCR). The amplified gene product is sequenced and the variation within the gene sequences is exploited to predict the taxonomic groups that are present in a sample, their relative abundances, and community-diversity measures.

These community descriptors are generated by a series of computational transformations of the original sequence data. Some of these transformations, such as sequence quality filtering, sequence alignments, and phylogeny building, are common bioinformatic tasks that can be accomplished with more general tools. Other transformations, such as taxonomic classification, or the

quantification of community-profile similarity are more likely to require databases and tools designed specifically for marker-gene analyses. Dozens of software tools and dependencies are required for a complete analysis, written in various programming languages and with documentation spread across many different locations. These tools interact with one another in a process workflow, feeding from one to the next. Often, the output format of one step does not match the input requirements of the next step, so an additional transformation is required. This can lead to complicated analyses that are performed with ad-hoc commands and data manipulations that make analysis replication and data-provenance tracking difficult or impossible.

This chapter demonstrates a microbial marker-gene analysis using the Quantitative Insights Into Microbial Ecology version 2 (QIIME2, pronounced “chime two”) software suite [3]. QIIME2 provides a software environment, data standards, and tool wrappers that allow for seamless interoperability between tools used for microbial community analysis. We describe a typical analysis pipeline using QIIME2, and demonstrate how study replication and data provenance can be simplified with scripting and QIIME artifacts. This chapter is accompanied by a GitHub repository (https://github.com/beiko-lab/mimb_16S) which contains scripts to download an example data set and process the data using the marker-gene analysis pipeline described here.

2 Materials

2.1 Sequence Data

This protocol requires a marker-gene data set generated from a 16S rRNA gene fragment and sequenced as paired-end reads on an Illumina platform. While in principle other genes and sequencing platforms could be used with QIIME2 and its associated tools, the default parameters and databases are tuned for the 16S rRNA gene and Illumina paired-end sequences. Use of sequence data from other genes or sequencing platforms would necessitate substituting an appropriate reference sequence set and a critical re-evaluation of default parameters and models on many steps, but particularly those associated with sequence denoising and taxonomic classification.

Sequence data should be in FASTQ format and must be named using the Illumina naming convention. For example, a gzip-compressed FASTQ file may be called `SampleName_S1_L001_R2_001.fastq.gz`, where `SampleName` is the name of the sample, `S1` indicates the sample number on the sample sheet, `L001` indicates the lane number, `R2` indicates that the file contains the reverse reads (with `R1` indicating forward reads), and the last three numbers are always `001` by convention. The files should be demultiplexed, which means that there is one FASTQ sequence file for every sample, and all of the FASTQ files should be placed

into the same directory. This directory should contain only two other files. The first is `metadata.yaml`. This is a simple text file that contains only the text `{phred-offset: 33}` on a single line (see **Note 1**). The second file is named `MANIFEST` and is a three-column comma-separated text file with the first column listing the sample name (matching the FASTQ filename convention), the second column listing the FASTQ file name, and the third column listing whether the reads are “forward” or “reverse.” Here are the first 5 lines of the `MANIFEST` file from the example data set:

```
sample-id,filename,direction
SRR3202913,SRR3202913_S0_L001_R1_001.fastq.gz,forward
SRR3202913,SRR3202913_S0_L001_R2_001.fastq.gz,reverse
SRR3202914,SRR3202914_S0_L001_R1_001.fastq.gz,forward
SRR3202914,SRR3202914_S0_L001_R2_001.fastq.gz,reverse
```

To import a directory of sequence files into QIIME, the directory must contain all of the FASTQ files, a `metadata.yaml` file, and a complete `MANIFEST` file listing each of the FASTQ files in the directory. Subheading 3.1 describes the import process.

2.2 Sample Metadata

Sample metadata is stored in a tab-separated text file. Each row represents a sample, and each column represents a metadata category. The first line is a header that contains the metadata category names. These cannot contain special characters and must be unique. The first column is used for sample names and must use the same names as in the `sample-id` column of the `MANIFEST` file. The QIIME developers host a browser-based metadata validation tool, Keemei (<https://keemei.qiime2.org/>), that checks for correct formatting and helps identify any errors in the metadata file [16]. Metadata files used in QIIME1 analyses are compatible with QIIME2 and can be used without modification.

2.3 Software

For this computational pipeline, we will be using the 2018.2 distribution of the QIIME2 software suite. The installation process has been significantly simplified over previous iterations of QIIME. The entire package, including all dependencies and tools, can be automatically installed with the Anaconda/Miniconda package and environment manager (available at <https://anaconda.org/>). The QIIME software is placed in a virtual environment so that it does not interfere or conflict with any existing software on the system. Once installed, the environment must be activated with the command `source activate qiime2-2018.2`, giving access to QIIME as well as the tools that it wraps. It is important to note that changes to the command line interface can occur between QIIME2 releases and with plugin updates. The companion GitHub repository will list any necessary changes to the protocol that may arise over time.

3 Methods

3.1 Import Data

In QIIME2, there are two main input/output file types: QIIME artifacts (.qza) and QIIME visualizations (.qzv). QIIME artifacts encapsulate the set of (potentially heterogeneous) data that results from a given step in the pipeline. The artifact also contains a variety of metadata including software versions, command parameters, timestamps, and run times. QIIME visualization files are analysis endpoints that contain the data to be visualized along with the code required to visualize it. Visualizations can be launched in a web browser with the `qiime tools view` command, and many feature interactive elements that facilitate data exploration. Data can be extracted from .qza or .qzv files using the `qiime tools extract` command (*see Note 2*).

If the sequences are in a directory named `sequence_data` (along with a `metadata.yml` and `MANIFEST` file, as described in Subheading 2.1), then the command to import these sequences into a QIIME artifact is:

```
qiime tools import --type
'SampleData[PairedEndSequencesWithQuality]' --input-path
sequence_data --output-path reads
```

The data type is specified as `SampleData[PairedEndSequencesWithQuality]`. This is QIIME's way of indicating that there are paired forward/reverse FASTQ sequence files for each sample (*see Note 3*). An output artifact named `reads.qza` will be created, and this file will contain a copy of each of the sequence data files (*see Note 4*).

3.2 Visualize Sequence Quality

The quality profile of sequences can vary depending on sequencing platform, chemistry, target gene, and many other experimental variables. The sequence qualities inform the choices for some of the sequence-processing parameters, such as the truncation parameters of the DADA2 denoising step [2]. QIIME includes an interactive sequence quality plot, available in the “q2-demux” plugin. The following command will sample 10,000 sequences at random and plot box plots of the qualities at each base position:

```
qiime demux summarize --p-n 10000 --i-data reads.qza
--o-visualization qual_viz
```

The plots, contained within a QIIME visualization .qzv file, can be viewed in a web browser by providing the .qzv file as an argument to the `qiime tools view` command. The sample size should be set sufficiently high to ensure an accurate representation of the qualities, but the run time of this command will increase with the sample size.

Since the `reads.qza` file was created from paired-end reads, the visualization will automatically display the quality distributions for a random sample of both the forward and reverse sequences. With Illumina paired-end data it is expected for there to be a decrease in the quality at the higher base positions. The point at which the quality begins to decrease should inform the truncation parameter used in the subsequent sequence denoising step. The truncation value is provided separately for forward and reverse sequence reads, so it is important to note where the quality decrease occurs for both sets.

3.3 Denoise Sequences With DADA2

As an alternative to OTU clustering at a defined sequence-identity cut-off (e.g., 97%), QIIME2 offers Illumina sequence denoising via DADA2 [2]. The `qiime dada2 denoise-paired` will both merge and denoise paired-end reads. The command has two required parameters: `--p-trunc-len-f` indicates the position at which the forward sequence will be truncated and `--p-trunc-len-r` indicates the position at which the reverse read will be truncated. Optional parameters include `--p-max-ee` which controls the maximum number of expected errors in a sequence before it is discarded (default is 2), and `--p-trunc-q` which truncates the sequence after the first position that has a quality score equal to or less than the provided value (default is 2). DADA2 requires the primers to be removed from the data to prevent false positive detection of chimeras as a result of degeneracy in the primers. If primers are present in the input sequence files, the optional `--p-trim-left-f` and `--p-trim-left-r` parameters can be set to the length of the primer sequences in order to remove them before denoising. The denoising process outputs two artifacts: a table file and a representative sequence file. The table file can be exported to the Biological Observation Matrix (BIOM) file format (an HDF5-based standard) using the `qiime tools export` command for use in other utilities [11]. The representative sequence file contains the denoised sequences, while the table file maps each of the sequences onto their denoised parent sequence.

3.4 Filter Sequence Table

After denoising with DADA2, many reads may have been excluded because they could not be merged or were rejected during chimera detection. You may wish to exclude any samples that have significantly fewer sequences than the majority. The `qiime feature-table summarize` command produces a visualization file that shows the spread of sequence depths across the samples. Use this visualization to identify a lower bound on the sequence depth and (if desired) filter out low sequence depth samples with the `qiime feature-table filter-samples` command with the `--p-min-frequency` parameter.

3.5 *Taxonomic Classification*

The QIIME2 software leverages the machine learning Python library scikit-learn to classify sequences [14]. A reference set can be used to train a naïve Bayes classifier which can be saved as a QIIME2 artifact for later re-use. This avoids re-training the classifier between runs, decreasing the overall run time. The QIIME2 project provides a pre-trained naïve Bayes classifier artifact trained against Greengenes (13_8 revision) trimmed to contain only the V4 hypervariable region and pre-clustered at 99% sequence identity [12]. To train a naïve Bayes classifier on a different set of reference sequences, use the `qiime feature-classifier fit-classifier-naive-bayes` command. Other pre-trained artifacts are available on the QIIME2 website (<https://docs.qiime2.org/>). Once an appropriate classifier artifact has been created or obtained, use the `qiime feature-classifier classify` command to generate the classification results.

3.6 *Visualize Taxonomic Classifications*

The taxonomic profiles of each sample can be visualized using the `qiime taxa barplot` command. This generates an interactive bar plot of the taxa present in the samples, as determined by the taxonomic classification algorithm and reference sequence set used earlier. Bars can be aggregated at the desired taxonomic level and sorted by abundance of a specific taxonomic group or by metadata groupings. Color schemes can also be changed interactively, and plots and legends can be saved in vector graphic format.

3.7 *Build Phylogeny*

A phylogenetic tree must be created in order to generate phylogenetic diversity measures such as unweighted and weighted UniFrac [9, 10] or Faith's phylogenetic diversity (PD) [7]. The process is split into four steps: multiple sequence alignment, masking, tree building, and rooting. QIIME2 uses MAFFT for the multiple sequence alignment via the `qiime alignment mafft` command [8]. The masking stage will remove alignment positions that do not contain enough conservation to provide meaningful information (default 40%) and can also be set to remove positions that are mostly gaps. The `qiime alignment mask` command provides this functionality. The tree building stage relies on FastTree (see **Note 5**) and can be invoked with `qiime phylogeny fasttree` [15]. The final step, rooting, takes the unrooted tree output by FastTree and roots it at the midpoint of the two leaves that are the furthest from one another. This is done using the `qiime phylogeny midpoint-root` command. The end result is a rooted tree artifact file that can be used as input to generate phylogenetic diversity measures.

3.8 *Compute Diversity Measures*

An array of alpha- and beta-diversity measures can be generated with a single command with QIIME2. The `qiime diversity core-metrics-phylogenetic` command will produce both phylogenetic and non-phylogenetic diversity measures, as well as alpha- and beta-diversity measures. As input, this command

requires a sequence/OTU table, a phylogenetic tree, and a sampling depth for random subsampling. A good value for the sampling depth is the number of sequences contained in the sample with the fewest sequences. It can be found by visualizing the `table_summary_output.qzv` file from the `qiime feature-table summarize` command. The `qiime diversity core-metrics-phylogenetic` command generates Faith's phylogenetic diversity, Shannon diversity, evenness, and observed OTUs (see **Note 6**) as alpha-diversity measures and weighted/unweighted UniFrac, Bray-Curtis, and Jaccard as beta-diversity measures. For each of the beta-diversity measures, QIIME2 automatically generates principal co-ordinates analysis visualizations. These are three-dimensional visualizations of the high-dimensional pairwise distance (or dissimilarity) matrices. These plots allow the researcher to identify groupings of similar samples at a glance.

3.9 Test for Diversity Differences Between Groups

We can test for significant differences between different sample groups using the `qiime diversity alpha-group-significance` and `qiime diversity beta-group-significance` commands. The alpha-diversity group significance command creates boxplots of the alpha-diversity values and significant differences between groups are assessed with the Kruskal-Wallis test. The beta-diversity command uses boxplots to visualize the distance between samples aggregated by groups specified in the metadata table file. Significant differences are assessed using a PERMANOVA analysis [1] or optionally with ANOSIM [5].

3.10 Alpha Rarefaction

An alpha rarefaction analysis is used to determine if an environment has been sequenced to a sufficient depth. This is done by randomly subsampling the data at a series of sequence depths and plotting the alpha diversity measures computed from the random subsamples as a function of the sequencing depth. A plateau on the rarefaction curve of a given sample provides evidence that the sample has been sequenced to a sufficient depth to capture the majority of taxa. Use the `qiime diversity alpha-rarefaction` command to generate the visualization file.

3.11 Data Provenance

Experimenting with different parameters and plugins can result in an accumulation of output files. It can be quite easy to lose track of which commands and parameters were used for which file. Thankfully, QIIME2 tracks the provenance of each artifact and visualization file. Using the online viewer (<https://view.qiime2.org>), an artifact or visualization file can be imported, presenting a provenance tab that shows the history of the file. The viewer lists each of the parameters used to create the file as well as the run time for the command and a comprehensive list of plugin and software versions. This information is provided not only for the imported file, but also for each of the files that were provided as input, and the input to

those files, and back until the data import command. This means that each QIIME artifact comes bundled with all the knowledge of how it was created.

4 Example

4.1 Retrieve Example Sequence Data

For the example analysis, we will be retrieving data from a recent study on the gut microbiome of the bumblebee, *Bombus pascuorum* [13]. In this study, 106 samples were collected from four different types of bumblebees. Twenty-four samples were collected from larvae (La), 47 from nesting bees (Nu), 18 from foragers that lived in the nest (Fn), 16 from foragers that were collected from the nearby environment (Fo), and one from the queen (Qu). The DNA from the microbiota in the midgut and hindgut of each insect was extracted and amplified using the 515f/806r primer pair (16S rRNA gene V4 region). DNA sequencing was performed on an Illumina MiSeq, resulting in a set of paired-end 16S rRNA gene fragment sequences with an insert size of approximately 254 bp in length.

The data were deposited in the European Bioinformatics Institute Short Read Archive (EBI SRA) at project accession PRJNA313530. The GitHub repository that accompanies this chapter (https://github.com/beiko-lab/mimb_16S) contains a BASH script named `fetchFastq.sh`. This script automatically downloads the raw FASTQ files and sample metadata from the EBI SRA.

4.2 Import Data

The `fetchFastq.sh` script creates a directory named `sequence_data/import_to_qiime/` that contains the forward and reverse FASTQ data files, the `MANIFEST` file, and the `metadata.yml` file. While in the directory containing the `fetchFastq.sh` script, the following command will import the FASTQ files into a QIIME artifact named `reads.qza`:

```
qiime tools import --type
'SampleData[PairedEndSequencesWithQuality]' --input-path
sequence_data/import_to_qiime --output-path reads
```

4.3 Visualize Sequence Quality

To generate visualizations of the sequence qualities, we run the command:

```
qiime demux summarize --p-n 10000 --i-data reads.qza
--o-visualization qual_viz
```

Next, the command `qiime tools view qual_viz.qzv` will launch the visualization in a web browser. Figure 1 shows the quality profile across a sample of 10,000 reverse reads. The quality




Fig. 1 Quality score box plots sampled from 10,000 random reverse reads

scores begin slightly lower, which is expected from the bases that belong to the primer sequences (*see* **Note 7**). We will be removing the primer sequences at the denoising stage. While the median quality scores remain fairly stable, the variance of the quality scores increases around position 130. We will trim some of the bases at the 5' end of the reverse reads during the denoising stage.

4.4 Denoise Sequences With DADA2

The example sequence data are 2×151 bp paired-end reads from an Illumina MiSeq using the 515f/806r primer set [4]. The quality plots (e.g., Fig. 1) indicate that the primers should be trimmed. The forward primer is 19 bp in length and the reverse primer is 20 bp in length, informing our choice for the parameters `--p-trim-left-f` and `--p-trim-left-r`, respectively. We see an increase in the variance of the quality scores for the reverse reads, so we will truncate the reverse reads at position 140. We will not truncate the forward reads, as the same dramatic increase in quality variance is not observed. Therefore, the `--p-trunc-len-f` parameter will be set to 151, and the `--p-trunc-len-r` parameter will be set to 140. At an average amplicon length of 254 bp, trimming the reverse reads by 11 bp would leave an average of 37 bp of overlap. This is sufficient for DADA2, which requires a minimum of 20 bp of overlap for the read merging step. Using `--p-n-threads 4` allows the program to perform parallel computations on 4 threads, and the `--verbose` option displays the DADA2 progress in the terminal.

```
qiime dada2 denoise-paired --i-demultiplexed-seqs reads.qza
--o-table table --o-representative-sequences
representative_sequences --p-trunc-len-f 151 --p-trunc-len-r 140
--p-trim-left-f 19 --p-trim-left-r 20 --p-n-threads 4 --verbose
```

The information printed to the terminal with the `--verbose` option shows a sampling of the sequence counts at each stage of the denoising process. In order to view the number of successfully denoised sequences for each sample, we create a summary of the output table file:

```
qiime feature-table summarize --i-table table.qza
--o-visualization table_summary
```

The visualization file provides detailed information about the denoised sequence counts, including the number of sequences per sample. Sample SRR3203007 had the fewest sequences, with 3615 non-chimeric denoised sequences identified by DADA2. The sample with the second-lowest sequencing depth was SRR3203003 with 42,138 sequences.

4.5 *Filter Sequence Table*

Since SRR3203007 has a significantly lower sequencing depth than all of the other samples, we will remove it from the table and exclude it from further analysis.

```
qiime feature-table filter-samples --i-table table.qza
--p-min-frequency 5000 --o-filtered-table filtered_table
```

This removes samples with fewer than 5000 sequences, which will remove only sample SRR3203007.

4.6 *Taxonomic Classification*

First, we must download the trained naïve Bayes classifier artifact. We will fetch this from the QIIME website with the command `wget`:

```
wget
https://data.qiime2.org/2018.2/common/gg-13-8-99-515-806-nb-
classifier.qza
```

This classifier artifact is trained on Greengenes August 2013 revision, trimmed to the V4 hypervariable region using primers 515f/806r, and clustered at 99% sequence identity. We then instruct QIIME to classify using this classifier artifact and the `scikit-learn` Python library:

```
qiime feature-classifier classify-sklearn --i-classifier
gg-13-8-99-515-806-nb-classifier.qza --i-reads
representative_sequences.qza --o-classification taxonomy
```




Fig. 2 Taxonomic profiles for the bumblebee gut samples at phylum level. Samples are sorted first by bee type and then by the abundance of the phylum *Proteobacteria*. Bee types are larval samples (La), nesting bees (Nu), foragers from the nest (Fn), foragers from the environment (Fo), and the queen (Qu)

4.7 Visualize Taxonomic Classifications

QIIME can generate interactive bar plots of the taxonomic profiles. The profiles can be sorted by metadata categories, so we have to provide the tab-separated metadata file which was generated by the data download script. The metadata file for the example data is located in `sequence_data/METADATA.txt`.

```
qiime taxa barplot --i-table filtered_table.qza --i-taxonomy
taxonomy.qza --m-metadata-file sequence_data/METADATA.txt
--o-visualization taxa-bar-plots
```

The resulting `taxa-bar-plots.qzv` visualization file can be launched with the command `qiime tools view taxa-bar-plots.qzv`. Figure 2 shows taxonomic profiles sorted by `bee_type` and then by the relative abundance of the phylum *Proteobacteria*.

4.8 Build Phylogeny

We build a phylogeny based on the 16S rRNA gene fragments in the four-step process described in Subheading 3.7. First, we align the denoised sequences with MAFFT.

```
qiime alignment mafft --i-sequences representative_sequences.qza
--o-alignment aligned_representative_sequences
```

Next, we mask the uninformative positions.

```
qiime alignment mask --i-alignment
aligned_representative_sequences.qza --o-masked-alignment
masked_aligned_representative_sequences
```

We build the phylogeny using the FastTree method.

```
qiime phylogeny fasttree --i-alignment
masked_aligned_representative_sequences.qza --o-tree unrooted_
tree
```

Finally, we root the tree at the midpoint, producing the `rooted_tree.qza` artifact that will be used as input to generate phylogenetic-diversity measures.

```
qiime phylogeny midpoint-root --i-tree unrooted_tree.qza
--o-rooted-tree rooted_tree
```

4.9 Compute Diversity Measures

We can use a single command to generate a series of phylogenetic and non-phylogenetic diversity measures. In order to compare samples with uneven sequencing depth, QIIME2 randomly subsamples or “rarefies” the sequences present in each environmental sample, at a user-specified depth. After filtering, our sample with the fewest sequences is SRR3203003 with 42,138 sequences. The smallest number of sequences in a given sample can be used as the subsampling depth, but here we will go slightly lower and use a depth of 41,000.

```
qiime diversity core-metrics-phylogenetic --i-phylogeny
rooted_tree.qza --i-table filtered_table.qza --p-sampling-depth
41000 --output-dir diversity_41000 --m-metadata-file
sequence_data/METADATA.txt
```

This command will generate several QIIME artifact files that contain the Faith’s phylogenetic diversity, observed OTUs, Shannon diversity, and evenness alpha-diversity measures for each sample. It also generates beta-diversity distance matrices for the Bray-Curtis, Jaccard, unweighted UniFrac, and weighted UniFrac measures, as well as visualizations of the principal co-ordinates analyses based on these distance measures. These visualization files have the filename suffix `_emperor.qzv` and when viewed will display a three-dimensional ordination plot. A static example of the interactive principal co-ordinates visualization is shown in Fig. 3.

4.10 Test for Diversity Differences Between Groups

We will test for significant differences in the microbial community diversity measures of the bee types. QIIME2 will perform the statistical tests for each of the sample groupings present in the metadata file. To run the tests for Faith’s phylogenetic diversity, we run:

```
qiime diversity alpha-group-significance --i-alpha-diversity
diversity_41000/faith_pd_vector.qza --m-metadata-file
sequence_data/METADATA.txt --o-visualization
diversity_41000/alpha_PD_significance
```




Fig. 3 A two-dimensional principal co-ordinates analysis ordination of the bee gut samples based on unweighted UniFrac distances. Samples are colored by bee type (blue: foragers from the environment, red: foragers from the nest, green: nesting bees, purple: larvae, orange: queen). Legend is visible in the “Color” tab of the interactive visualization. The “Axis” tab, shown here, allows any of the principal co-ordinates to be selected, allows the background and font colors to be changed, and shows the Scree plot that indicates the proportion of the data variation captured by each principal co-ordinate

This runs all-group and pairwise Kruskal-Wallis tests, a non-parametric analysis of variance. The visualization file, `alpha_PD_significance.qzv`, presents boxplots and test statistics for each metadata grouping. In this example data set, the larval samples had a significantly higher phylogenetic diversity value compared with each of the other bee types ($p < 0.001$).

The test for beta-diversity group distances can be performed with PERMANOVA (the default method) or ANOSIM. The `qiime diversity beta-group-significance` command computes only one metadata grouping at a time, so to test the differences between bee types we have to supply the appropriate column name from the metadata file:

```
qiime diversity beta-group-significance --i-distance-matrix
diversity_41000/bray_curtis_distance_matrix.qza --m-metadata-file
sequence_data/METADATA.txt --m-metadata-column bee_type
--o-visualization diversity_41000/beta_bray_beetype_significance
```

This runs an all-group PERMANOVA analysis on the Bray-Curtis dissimilarity measures for each bee type. For this data set, the PERMANOVA test reveals that the five bee types have significant differences in their community compositions. The pairwise distance boxplots show that this is largely driven by the larval samples, an observation corroborated by the relatively robust clustering of larval samples visible in the principal co-ordinates ordination (Fig. 3).

4.11 Alpha Rarefaction

Our final analysis is to create alpha rarefaction curves in order to determine if the samples have been sequenced to a sufficient depth. The `qiime diversity alpha-rarefaction` command will generate rarefaction curves based on the Shannon diversity and observed OTUs measures by default, and will additionally generate phylogenetic diversity-based curves if the phylogenetic tree created above is provided using the `--i-phylogeny` parameter. The desired alpha-diversity measure is selected interactively after the visualization file is launched.

```
qiime diversity alpha-rarefaction --i-table filtered_table.qza
--p-max-depth 41000 --o-visualization
diversity_41000/alpha_rarefaction.qzv --m-metadata-file
sequence_data/METADATA.txt --i-phylogeny rooted_tree.qza
```

Figure 4 (*see Note 8*) shows the alpha rarefaction curves with the results average by bee type. We can derive a few insights from this table. The first is that each of the bee type categories appear to plateau. Although the diversity measure does generally continue to increase as a function of the sequencing depth, the accumulation slows significantly, suggesting that we have sufficient sample sequence depth to have captured the majority of taxa present in the sample. New taxa that may be picked up by additional sequencing effort are likely to be either rare microorganisms or the result of sequencing error. The second thing we can learn from Fig. 4 is that the larval samples have a significantly higher phylogenetic diversity than the other bee types. This result agrees with the Kruskal-Wallis test performed in Subheading 4.10.

4.12 Data Provenance

Marker-gene analyses have the potential to generate many output files. Data-provenance tracking ensures that we do not lose track of how each file was generated. If you forget or are unsure how a file




Fig. 4 Rarefaction curves based on the phylogenetic diversity measure, with samples aggregated by bee type

was generated, it can be imported to the viewer at <https://view.qiime2.org> and the history is available in the “Provenance” tab. For example, Fig. 5 shows the provenance graph for the unweighted UniFrac principal co-ordinates file generated in Subheading 4.9, `unweighted_unifrac_emperor.qzv`. Clicking on the circles in the graph reveals the output file type, format, and unique identifier. Clicking on the blocks reveals the command that generated the output file(s). The arrows show the flow of output files from one step as input to the next.




Fig. 5 Data provenance graph for the `unweighted_unifrac_emperor.qzv` file generated in the example analysis

5 Notes

1. Most recent Illumina sequence data is produced with the CASAVA pipeline 1.8+ which uses the PHRED+33 encoding indicated by the “metadata.yml” text `{phred-offset: 33}`. If your FASTQ data was generated by an older version of the CASAVA pipeline, the qualities will be in PHRED+64 format. In this case, the “metadata.yml” text should be changed to `{phred-offset: 64}`.
2. The .qza and .qvz files are simply re-named ZIP files. The file extension can be changed from .qza or .qvz to .zip and extracted with any ZIP file decompression tool on systems where QIIME is not installed or available, such as Windows PCs.
3. A list of the data types in QIIME2 (2018.2 distribution) is available at <https://docs.qiime2.org/2018.2/semantic-types/>. Single-end FASTQ files correspond to the data type “SampleData[SequencesWithQuality].”
4. Preserving copies of the raw input data within an artifact enhances research reproducibility by allowing the source data to be easily identified and ensures that the data are more closely coupled to the analysis. Just keep an eye on your hard drive usage!
5. The QIIME2 2018.2 distribution comes with FastTree version 2.1.10 compiled with double precision. This mitigates issues with resolving short branch lengths that could occur when using the version of FastTree that was distributed with earlier versions of the QIIME software suite.
6. Even though we are not using the common 97% OTU clustering approach, the denoised sequences from DADA2 can still be considered operational taxonomic units, so the name “observed OTUs” is still appropriate for this diversity measure. However, to avoid confusion, you may wish to describe the “observed OTUs” measure as “observed taxa” or “observed representative sequences.”
7. The lower quality scores at the beginning of each read are caused by the homogeneity of the primer sequences. This makes it difficult for the Illumina sequencer to properly identify clusters of DNA molecules [6].
8. Several of the QIIME visualizations do not have a clear way to export the plots in high-quality scalable vector graphic (SVG) format. The alpha rarefaction plot is one such visualization. A simple screenshot will result in a raster graphic that may have sufficient resolution for inclusion in a research article. The New York Times’ SVG Crowbar utility (<https://nytimes.>

github.io/svg-crowbar/) can be used to extract many of these plots in SVG format. These files can then be easily manipulated using Inkscape or similar vector graphic-editing programs.

References

1. Anderson M (2005) PERMANOVA: a FORTRAN computer program for permutational multivariate analysis of variance, 24th edn. Department of Statistics, University of Auckland, Auckland
2. Callahan BJ, McMurdie PJ, Rosen MJ, Han AW, Johnson AJA, Holmes SP (2016) DADA2: high-resolution sample inference from Illumina amplicon data. *Nat Methods* 13(7):581
3. Caporaso JG, Kuczynski J, Stombaugh J, Bittinger K, Bushman FD, Costello EK, Fierer N, Pena AG, Goodrich JK, Gordon JI et al (2010) QIIME allows analysis of high-throughput community sequencing data. *Nat Methods* 7(5):335
4. Caporaso JG, Lauber CL, Walters WA, Berg-Lyons D, Lozupone CA, Turnbaugh PJ, Fierer N, Knight R (2011) Global patterns of 16S rRNA diversity at a depth of millions of sequences per sample. *Proc Natl Acad Sci USA* 108:4516–4522
5. Chapman M, Underwood A (1999) Ecological patterns in multivariate assemblages: information and interpretation of negative values in ANOSIM tests. *Mar Ecol Prog Ser* 180:257–265
6. Fadrosch DW, Ma B, Gajer P, Sengamalai N, Ott S, Brotman RM, Ravel J (2014) An improved dual-indexing approach for multiplexed 16S rRNA gene sequencing on the Illumina MiSeq platform. *Microbiome* 2(1):6
7. Faith DP (1992) Conservation evaluation and phylogenetic diversity. *Biol Conserv* 61(1):1–10
8. Katoh K, Standley DM (2013) MAFFT multiple sequence alignment software version 7: improvements in performance and usability. *Mol Biol Evol* 30(4):772–780
9. Lozupone C, Knight R (2005) UniFrac: a new phylogenetic method for comparing microbial communities. *Appl Environ Microb* 71(12):8228–8235
10. Lozupone CA, Hamady M, Kelley ST, Knight R (2007) Quantitative and qualitative β diversity measures lead to different insights into factors that structure microbial communities. *Appl Environ Microb* 73(5):1576–1585
11. McDonald D, Clemente JC, Kuczynski J, Rideout JR, Stombaugh J, Wendel D, Wilke A, Huse S, Hufnagle J, Meyer F et al (2012) The Biological Observation Matrix (BIOM) format or: how I learned to stop worrying and love the ome-ome. *GigaScience* 1(1):7
12. McDonald D, Price MN, Goodrich J, Nawrocki EP, DeSantis TZ, Probst A, Andersen GL, Knight R, Hugenholtz P (2012) An improved Greengenes taxonomy with explicit ranks for ecological and evolutionary analyses of bacteria and archaea. *ISME J* 6(3):610
13. Parmentier A, Meeus I, Nieuwerburgh F, Deforce D, Vandamme P, Smagghe G (2018) A different gut microbial community between larvae and adults of a wild bumblebee nest (*Bombus pascuorum*). *Insect Sci* 25(1):66–74
14. Pedregosa F, Varoquaux G, Gramfort A, Michel V, Thirion B, Grisel O, Blondel M, Prettenhofer P, Weiss R, Dubourg V et al (2011) Scikit-learn: machine learning in Python. *J Mach Learn Res* 12(Oct):2825–2830
15. Price MN, Dehal PS, Arkin AP (2010) FastTree 2—approximately maximum-likelihood trees for large alignments. *PLoS One* 5(3):e9490
16. Rideout JR, Chase JH, Bolyen E, Ackermann G, González A, Knight R, Caporaso JG (2016) Keemei: cloud-based validation of tabular bioinformatics file formats in Google Sheets. *GigaScience* 5(1):27